


74. ***COLLAR STIFF NECK REGULAR***

Specifications

| Qualitative Requirements | Picture of Equipment |
|---|--|
| <p>COLLAR STIFF NECK REGULAR(AS PER IS: 11569)</p> <ul style="list-style-type: none">a) Rigid cervical collar-extrication typeb) Material: soft plastic rigid enough to avoid 70% of neck movement and lined inside with soft cushion.c) Padded chin rest for application by sweep technique.d) Colour coded and measure holes for easy identification and measurement of measured of correct size.e) Velcro for easy fixation.f) MRI/CT Scan compatible; radio lucentg) Tracheal access opening for advanced airway procedure;h) Rear panel opening for cervical spine palpationi) Latex free and durablej) This size of finished c-collar should not be less than 6 cm when measurement as per prescribed procedure. |  |