

77. **COLLAR STIFF NECK TALL**

Specifications

Qualitative Requirements	Picture of Equipment
<p>(AS PER IS: 11569)</p> <ul style="list-style-type: none">a) Rigid cervical collar- extrication typeb) Material: soft plastic rigid enough to avoid 70% of neck movement and lined inside with soft cushion.c) Padded chin rest for application by sweep technique.d) Colour coded and measure holes for easy identification and measured of correct size.e) Velcro for easy fixation.f) MRI/CT Scan compatible; radiolucentg) Tracheal access opening for advanced airway procedure;h) Rear panel opening for cervical spine palpationi) Latex free and durable.j) This size finished c-collar should not be less than 10 cm when measurement as per prescribed procedure.	