

75. **COLLAR STIFF NECK ' NO NECK'**

Specifications

Qualitative Requirements

COLLAR STIFF NECK (NO NECK) (AS PER IS: 11569)

- a) Rigid cervical collar- extrication type
- b) Material: soft plastic rigid enough to avoid 70% of neck movement and lined inside with soft cushion.
- c) Padded chain rest for application by sweep technique.
- d) Colour coded and measure holes for easy identification and measurement of correct size.
- e) Velcro for easy fixation.
- f) MRI/CT Scan compatible; radio lucent
- g) Tracheal access opening for advanced airway procedure;
- h) Rear panel opening for cervical spine palpation
- g) Latex free and durable
- i) The size of finished c-collar should not be less than 3 cm when measurement as per prescribed procedure.

Picture of Equipment

